

AL03-02-230 Thermostatic water bath 4.5 litre



Strength	Class III Thermostatic water bath +20 to +150C,4.5 litre, at attractive price
Size (liter)	4.5
Temperature C	20/150C
Control Panel	Digital
Bath Volume max. Litre	4.5
Bath Work Area (WxL/D)(mm)	150*150/150
Heating Perf. (KW)	2.0 kW
Pump Pressure bar	0.3 at 0l
Pump Flow L/min.	14l/min at 0bar
Voltage(V)(~10%) 50/60Hz	230
Width (mm)	170
Length (mm)	330
Height (mm)	360
Color	Cream white
Weight (empty)(kg)	8