

AL03-03-230 Thermostatic water bath 12 litre



Strength	Class III Thermostatic water bath 20 to 150C, 12 litre, at attractive price
Size (liter)	12
Temperature C	20/150C
Control Panel	Digital
Bath Volume max. Litre	12
Bath Work Area (WxL/D)(mm)	180*130/150
Heating Perf. (KW)	2.0 kW
Pump Pressure bar	0.3 at 0l
Pump Flow L/min.	14l/min at 0bar
Voltage(V)(~10%) 50/60Hz	230
Width (mm)	390
Length (mm)	330
Height (mm)	370
Color	Cream white
Weight (empty)(kg)	9