

Advantage-Lab

The best priced laboratory equipment

PRODUCT DESCRIPTION

AL06-02-230 UNIVERSAL CENTRIFUGE

Centrifuges to meet your needs

Maximum Volume 400ml (4x100ml)

Many rotors and adapters available for both swing-out and fixed angle rotors as well as micro plate rotors (see AL06-02-400/420/430/440 and adapted buckets)

Rotors are easy to exchange

Runs can be limited from 10 seconds up to 99 hours 59 seconds (in steps of 10 seconds up to 59 minutes and then steps of 1 minute) or can be on permanent run.

Encapsulated motor

Highly versatile centrifuge to fit your every need

Centrifuge Safety features

Rotor type recognition: manual input for adjustment before start up

Imbalance sensor with safety shutdown mechanism

Electronic Lid locking system prevents opening of lid while in usage

Error mention with imbalance or other failure

Manufactured in accordance to IEC 1010 regulations

Corrosion resistant steel chambers

Highly durable construction provides stability during imbalance

User friendly Operating features

Large easy to view control panel with digital display

Easy to adapt – first offers last used program

Control buttons for selecting

Speed /rcf

Running time

Acceleration – Deceleration

Digital displays for pre-set or actual values

Selectable acceleration and deceleration rates (braking can be adjusted to 10 different intensities – position 0 means no braking)